MASCARILLA QUIRURG TI 3/C 50UD

30/05/2020 V1.0

Código 018188

Descriptiva

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Specifications of Single-use Face Mask

Prepared by: Li Minglu/Sun Yiyuan

Checked by: Shen Haigen

Reviewed by: Guan Yubo

1 Introduction

The mask consists of three layers of non-woven fabric. The outer layer is blue polypropylene spunbond non-woven fabric, the middle layer is polypropylene meltblown non-woven fabric with bacteria and particle filtering effect, and the inner layer is white polypropylene spunbond non-woven fabric. The masks are suitable for wearing in an ordinary medical environment, blocking the exhalation or spraying of pollutants from the oral or nasal cavities. Used in non-invasive operation, bacterial filtration efficiency is over 95%.

2 Standard

The mask products are based on standards: EN 14683: 2005 Type I (CE)

3 Product appearance

Fig.1 product appearance

Fig.2 product appearance

4 Specification

4.1 Size and appearance requirements


Fig.3 mask structure

1. Mask length >170 mm
2. Mask width: 165 mm, three stacks in the middle, the width of the mask after stacked is 95 mm, and the width of each stack ≥ 10 mm. The nose-clip side as the upper part and the dark colour as the outer side. After opening the stack, the opening side should face down
3. The width of the pressing edge ≤ 10 mm, and the upper nose-clip ≤ 16 mm
4. Nose-clip: Length >80 mm width is preferably 3mm
5. Ear-loop: The length is 180 mm and diameter is 3 mm. The material is polyester spandex. The inner layer is welded and the distance from the mask edge is not more than 10 mm.

Table 1 Parameters of single-use face masks

Items	Parameters
Product name	Single-use face mask
Textures	Polypropylene spunbond non-woven fabric, polypropylene melt-blown non-woven fabric, metal core plastic nose-clip.
Model	Flat
Size	175 mm * 95 mm
Application	Used mainly for respiratory protection of related personnel, which can filter particles in the air and block droplets and microorganisms.
Expiration date	The validity period is subject to the validity period promised by the non-woven fabric supplier, and the maximum period shall not exceed three years.
Packaging size	5 pcs / bag, 50 pcs / box, 2000 pcs / carton

Items	Parameters
Sterilization type	Non-sterile
Usage	1. Flatten the mask and hang both sides of the mask on your ears; 2. Fix the nose-clip in the shape of the nose to prevent dirty air from entering; 3. Pull the mask to the lower jaw so as to produce a good protective effect.
Special storage conditions and methods	Store in a relative humidity not exceeding 80%, non-corrosive gas, in a well-ventilated room and avoiding high temperatures.
Remarks	Epidemic emergency products

4.2 Inner and outer spunbond non-woven fabrics

The inner and outer layers are polypropylene spunbond non-woven fabrics. The outer layer is generally blue, the inner layer is white and the areal density is 25g / m². The width is based on the equipment width, the inner layer width is 200 mm, and the outer layer width is 180 mm. The packaging method is coil and the diameter of the coil cannot be greater than the equipment mounting limit. It requires sealed packaging.

The technical requirements for the inner and outer spunbond nonwovens are shown in Table 2:

Table 2 the technical requirements for the inner and outer spunbond nonwovens

Items		Index
Areal density		25 g/m ²
Breaking strength	Longitudinal	26 N/5cm
	Horizontal	18 N/5cm
Fluorescence		None
Formaldehyde content		≤ 300 mg/kg
pH		5.4-7.6
Microbial index	Total number of bacteria	cfu/g ≤20
	Coliform	None
	Pathogenic pyogenic bacteria	None
	Total fungal colonies	None
Biological evaluation	Cytotoxicity	Not more than 2 levels
	Skin irritation	Primary stimulation index does not exceed 0.4
	Delayed hypersensitivity	None

4.3 Interlayer meltblown non-woven fabric

The meltblown non-woven fabric is made of high-speed hot air to draw polypropylene (PP) to form ultra-fine fiber agglomeration. Its specifications are white and its areal

density is 25 g / m². The packaging method is coil, the diameter of the coil cannot be larger than the equipment mounting limit and the package is sealed.

The technical requirements for the interlayer meltblown non-woven fabric are shown in Table 3

Table 3 the technical requirements for the interlayer meltblown non-woven fabric

Items		Index
Areal density		25 g/m ²
Fluorescence		None
Formaldehyde content		≤ 300 mg/kg
pH		5.4-7.6
Microbial index	Total number of bacteria	cfu/g ≤20
	Coliform	None
	Pathogenic pyogenic bacteria	None
	Total fungal colonies	None
Particle filtration efficiency		≥ 30 %
Bacterial filtration efficiency		≥ 95 %

4.4 Left and right reinforcing edges

The left and right reinforcing edges are polypropylene spunbond non-woven fabrics, which are used to fix the ear-loop and three layers of non-woven fabrics together. 20 mm wide, white, areal density 35 g / m², packed in rolls, sealed.

4.5 Nose-clip

Nose-clip is made of PP plastic with galvanized iron wire core, length ≥ 80 mm, width 2.5 mm ~ 3 mm. It has plasticity and fits the nose shape by pressing down to increase the effective tightness of the mask.

4.6 Ear loop

The ear loop is generally made of polyester spandex which is connected to the mask body by ultrasonic welding. The ear-loop elasticity is designed to fit comfortably and to effectively fix the position of the mask in place. Generally, it is 180 mm. The package is sealed.

5 Product testing requirements

According to Article (8) of the ‘Guiding Principles for the Technical Review of Medical Mask Product Registration’:

Product testing requirements: The testing of medical mask products includes factory inspection and type inspection. The factory inspection items should have at least the following items: appearance, structure and size, nose-clips, mask bands,

microbiological indicators, and ethylene oxide residue (if ethylene oxide sterilization is used). The type inspection shall be the full performance inspection of the product standard. Specific testing items and requirements are shown in Table 4.

Table 4 Testing items and requirements for finished masks

S.N.	Test items	Requirements
1	Appearance	Type inspection, factory inspection
2	Structure and size	Type inspection, factory inspection
3	Nose-clip	Type inspection, factory inspection
4	Ear loop	Type inspection, factory inspection
5	Ventilation resistance	Type inspection
6	Bacterial filtration efficiency (BFE)	Type inspection
8	Sterilization requirements	E.O. sterilization, irradiation sterilization is prohibited, non-sterile products are not required
9	Skin irritation	Type inspection, supplier raw material inspection report
10	Cytotoxicity	Type inspection, supplier raw material inspection report
11	Delayed hypersensitivity	Type inspection, supplier raw material inspection report

6 Requirements for production environment and process hygiene

The production workshop should meet the requirements of '100,000 air cleanliness class'.

Requirements:

1. Maximum allowable number of dust particles: no more than 3.5 million particles $\geq 0.5 \mu\text{m}$; no more than 20,000 particles $\geq 5 \mu\text{m}$;
2. The maximum allowable number of microorganisms: the number of planktonic bacteria does not exceed 500 / m^3 ; the number of sedimentary bacteria does not exceed 10 / culture dish;
3. Pressure difference: The pressure difference between clean rooms of the same cleanliness level is the same. The pressure difference between adjacent clean rooms of different cleanliness levels should be $\geq 5 \text{ Pa}$, and the clean room and non-clean rooms should be $\geq 10 \text{ Pa}$ (mainly to ensure that air from the clean area flow to non-clean areas to avoid airflow backwashing);
4. Temperature and humidity. The temperature is generally controlled $20 \sim 22 \text{ }^\circ\text{C}$ in winter; $24 \sim 26^\circ\text{C}$ in summer; fluctuation $\pm 2 \text{ }^\circ\text{C}$. The humidity of clean room in winter is controlled at 30-50%, and the humidity of clean room which in summer

is controlled at 50-70%. When there are no special requirements for temperature and humidity, it is advisable to wear clean work clothes without discomfort.

7 Labeling and packaging requirements

According to regulatory requirements, product packaging should include the following information.

7.1 Labeling requirements

1. Product name, size and model type;
2. Production date / lot number, manufacture date or expiration date;
3. Manufacturer name, address and contact information;
4. Name, address, and contact information of the registrant or recorder;
5. Implementation of standard number, product registration certificate number, production license number / production filing certificate number;
6. The record number of the equipment urgently needed for the epidemic situation;
7. Instructions for use (at least including front and back identification and wearing methods);
8. Storage conditions, special storage, operating conditions or instructions;
9. 'Single use' words or symbols;
10. If it is a sterilised product, it should have a corresponding sterilisation mark, and indicate the sterilisation method used and the sterilisation period
11. Product Usage;
12. The 'emergency product' is marked clearly;
13. If the medical device can be marketed as sterile or non-sterile, the packaging or labelling of the product should be distinguished;
14. According to product characteristics, graphics, symbols and other related content should be marked;
15. Necessary warnings, precautions;

7.2 Specification requirements:

1. Product name, model, specifications;
2. Applicant / manufacturing company name, residence, production address, contact information and after-sales service unit;
3. Recording certificate number of urgently needed equipment for prevention and control;
4. Product performance, main structural composition or composition, scope of application;
5. Contraindications, precautions, warnings and tips;

-
6. Installation and use instructions or diagrams, medical devices used by consumers themselves should also have special instructions for safe use;
 7. Product maintenance and maintenance methods, special storage and transportation conditions and methods;
 8. Production date or expiration date;
 9. List of spare parts, including parts, accessories, replacement cycle of consumables, and instructions for replacement methods;
 10. Explanation of graphics, symbols, abbreviations, etc. used in medical device labels;
 11. Date of preparation or revision of the manual;
 12. The 'emergency product' must be marked prominently.

7.3 Packaging requirements

The PE bag should be sealed and packed in a carton. Use as soon as possible after the sealed bag is opened to prevent contamination.

BYD Single-use Face Mask User Instructions

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
EN 14683:2005 Type I

Product Model: FLAT

BYD Single-use face masks consist of three layers of nonwoven material

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Middle Layer

Polypropylene melt-blown nonwoven with pathogen filtering

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
Nose clip
Barrette nasale
Clip per naso
Nasclip
Neusclip
Clip para la nariz

Product Specifications

- 1 Mask length: > 170 mm.
- 2 Mask expansion width: 165 mm. Post-stack width is 95 mm with three stacks in the middle. Each stack width is not less than 10 mm. The nose clip is located on the upper edge of the mask, and the outside of the mask is darker in color.
- 3 The width of the support part is not more than 10 mm, and the position of the upper nose clip is not more than 16 mm.
- 4 Metal core plastic nose-clip: length > 80 mm; width approximately 3 mm.
- 5 Ear loop: 180 mm in length and 3 mm in diameter, composed of polyester, spandex and other materials.

Parameter

Product name	Single-use Face Mask
Bacterial filtration efficiency	≥ 95 %
Model	Flat
Size	175mm x 95mm
Weight	3.2 g
Expiration date	2 years after production
Packaging specification	10 pcs/bag; 50 pcs/box; 2000 pcs/carton
Packaging dimensions (L x W x H)	Carton: 520 x 400 x 470mm Box: 190x100x110mm Bag: 220 x 108mm
Sterilization	Non-sterile

Important limitations

- 1 Please use the product within the period of validity.
- 2 One-time use only. Disposal of contaminated products should be considered as hazardous waste and in accordance with national regulations.
- 3 Duration of use: 4 hours.
- 4 Do not use if packaging is damaged.
- 5 Indesirable effects and incidations: in contact with the skin, the mask can cause allergic reactions on sensitive people. If this is the case, remove the mask and seek medical advice.
- 6 Please read the instructions carefully before using.

Storage

Store in a well-ventilated place with relative humidity below 80%; avoid high temperatures, corrosive gas and exposure to flame.

Warning

This mask does not eliminate the risk of contracting any disease or infection.
Change immediately if contaminated with blood or body fluid.

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
Wellkang Ltd
16 Castle St, Dover, Kent, CT16 1PW, England, UK

Tel: +86 755-89888888

The following sample(s) was/were submitted and identified on behalf of the applicant as:

Applicant: BYD Precision Manufacture Co., Ltd.
NO.3001 Baohe Road, Baolong Industrial Area, Longgang, Shenzhen

Manufacturer: The same as the applicant

Sample Name: Single-use face mask

Lot No. : /

Quantity: 50 PCS

Model/Type Reference: Flat/17.5x9.5cm

Date of Sample Received: 2020-03-19

Date of Testing: 2020-03-19 to 2020-04-06

The Standards: EN 14683:2019

Test Items: Bacterial filtration efficiency (BFE), Differential pressure and Microbial cleanliness as required by applicant

Test Result: -Please refer to next page(s)-

Remark: The test was carried out by external laboratory which assessed as competent.

Jova Lan

Jova Lan
Technical Manager
Medical Laboratory

SGS-CSTC Standards Technical Services Co., Ltd. Guangzhou Branch
Table 1

Unless otherwise agreed in writing, this document is issued by the Company subject to its General Conditions of Service printed overleaf, available on request or accessible at <http://www.sgs.com/en/Terms-and-Conditions.aspx> and, for electronic format documents, subject to Terms and Conditions for Electronic Documents at <http://www.sgs.com/en/Terms-and-Conditions/Terms-e-Document.aspx>. Attention is drawn to the limitation of liability, indemnification and jurisdiction issues defined therein. Any holder of this document is advised that information contained hereon reflects the Company's findings at the time of its intervention only and within the limits of Client's instructions, if any. The Company's sole responsibility is to its Client and this document does not exonerate parties to a transaction from exercising all their rights and obligations under the transaction documents. This document cannot be reproduced except in full, without prior written approval of the Company. Any unauthorized alteration, forgery or falsification of the content or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law. Unless otherwise stated the results shown in this test report refer only to the sample(s) tested and such sample(s) are retained for 30 days only.
Attention: To check the authenticity of testing / inspection report & certificate, please contact us at telephone: (86-755) 8307 1443, or email: CN.Doccheck@sgs.com
No.198 Kexu Road, Sciencetech Park, Guangzhou Economic & Technology Development District, Guangzhou, China 510663 t (86-20) 82155555 f (86-20) 82075058 www.sgsgroup.com.cn
中国·广州·经济技术开发区科学城科珠路198号 邮编: 510663 t (86-20) 82155555 f (86-20) 82075058 sgs.china@sgs.com

Test Items		Result - Remarks	Test Methods
Differential Pressure (Pa/cm ²)	1	24.8	EN 14683:2019 Annex C
	2	24.4	
	3	26.5	
	4	25.2	
	5	26.0	

Table 2

Test Items		Result - Remarks	Test Methods
Bacterial Filtration Efficiency (BFE)(%) Staphylococcus aureus ATCC 6538	1	99.1	EN 14683:2019 Annex B
	2	99.5	
	3	99.5	
	4	99.2	
	5	99.1	

Note: Control average: 2196 CFU.

Note: Mean particle size: 2.8 μ m.

Note: Testing side: outside of specimen.

Note: Testing area: 39.5 cm².

Unless otherwise agreed in writing, this document is issued by the Company subject to its General Conditions of Service printed overleaf, available on request or accessible at <http://www.sgs.com/en/Terms-and-Conditions.aspx> and, for electronic format documents, subject to Terms and Conditions for Electronic Documents at <http://www.sgs.com/en/Terms-and-Conditions/Terms-e-Documents.aspx>. Attention is drawn to the limitation of liability, indemnification and jurisdiction issues defined therein. Any holder of this document is advised that information contained hereon reflects the Company's findings at the time of its intervention only and within the limits of Client's instructions, if any. The Company's sole responsibility is to its Client and this document does not exonerate parties to a transaction from exercising all their rights and obligations under the transaction documents. This document cannot be reproduced except in full, without prior written approval of the Company. Any unauthorized alteration, forgery or falsification of the content or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law. Unless otherwise stated the results shown in this test report refer only to the sample(s) tested and such sample(s) are retained for 30 days only.
Attention: To check the authenticity of testing / inspection report & certificate, please contact us at telephone: (86-755) 8307 1443, or email: CN.Doccheck@sgs.com
No.198 Kezhu Road, Scientech Park, Guangzhou Economic & Technology Development District, Guangzhou, China 510663 t (86-20) 82155555 f (86-20) 82075058 www.sgsgroup.com.cn
中国·广州·经济技术开发区科学城科珠路198号 邮编: 510663 t (86-20) 82155555 f (86-20) 82075058 sgs.china@sgs.com

Table 3

Test Items		Result - Remarks	Test Methods
Microbial Cleanliness (cfu/g)	1	8.2	EN 14683:2019 EN ISO 11737-1:2018
	2	18.2	
	3	24.0	
	4	4.6	
	5	13.8	

Unless otherwise agreed in writing, this document is issued by the Company subject to its General Conditions of Service printed overleaf, available on request or accessible at <http://www.sgs.com/en/Terms-and-Conditions.aspx> and, for electronic format documents, subject to Terms and Conditions for Electronic Documents at <http://www.sgs.com/en/Terms-and-Conditions/Terms-e-Documents.aspx>. Attention is drawn to the limitation of liability, indemnification and jurisdiction issues defined therein. Any holder of this document is advised that information contained hereon reflects the Company's findings at the time of its intervention only and within the limits of Client's instructions, if any. The Company's sole responsibility is to its Client and this document does not exonerate parties to a transaction from exercising all their rights and obligations under the transaction documents. This document cannot be reproduced except in full, without prior written approval of the Company. Any unauthorized alteration, forgery or falsification of the content or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law. Unless otherwise stated the results shown in this test report refer only to the sample(s) tested and such sample(s) are retained for 30 days only.

Attention: To check the authenticity of testing / inspection report & certificate, please contact us at telephone: (86-755) 8307 1443, or email: CN.Doccheck@sgs.com

SGS-CTC Standards Technical Services Co., Ltd.
Guangzhou Branch / Testing Center / EEC Laboratory

No.198 Kezhu Road, Sciencetech Park, Guangzhou Economic & Technology Development District, Guangzhou, China 510663 t (86-20) 82155555 f (86-20) 82075058 www.sgsgroup.com.cn
中国·广州·经济技术开发区科学城科珠路198号 邮编: 510663 t (86-20) 82155555 f (86-20) 82075058 sgs.china@sgs.com

Photo documentation:

Sample description

Single-use face mask

-- End of this report--

SGS-CTC Standards Technical Services Co., Ltd.
Guangzhou Branch of Testing Center EEC Laboratory

Unless otherwise agreed in writing, this document is issued by the Company subject to its General Conditions of Service printed overleaf, available on request or accessible at <http://www.sgs.com/en/Terms-and-Conditions.aspx> and, for electronic format documents, subject to Terms and Conditions for Electronic Documents at <http://www.sgs.com/en/Terms-and-Conditions/Terms-e-Documents.aspx>. Attention is drawn to the limitation of liability, indemnification and jurisdiction issues defined therein. Any holder of this document is advised that information contained hereon reflects the Company's findings at the time of its intervention only and within the limits of Client's instructions, if any. The Company's sole responsibility is to its Client and this document does not exonerate parties to a transaction from exercising all their rights and obligations under the transaction documents. This document cannot be reproduced except in full, without prior written approval of the Company. Any unauthorized alteration, forgery or falsification of the content or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law. Unless otherwise stated the results shown in this test report refer only to the sample(s) tested and such sample(s) are retained for 30 days only.
Attention: To check the authenticity of testing / inspection report & certificate, please contact us at telephone: (86-755) 8307 1443, or email: CN.Doccheck@sgs.com

No.198 Kezhu Road, Sciencetech Park, Guangzhou Economic & Technology Development District, Guangzhou, China 510663 t (86-20) 82155555 f (86-20) 82075058 www.sgsgroup.com.cn
中国·广州·经济技术开发区科学城科珠路198号 邮编: 510663 t (86-20) 82155555 f (86-20) 82075058 sgs.china@sgs.com

Member of the SGS Group (SGS SA)

EC Declaration of Conformity

according to the Medical Devices Directive 93/42/EEC

Class I Medical Device

(non-sterile, without measuring function)

Manufacturer: BYD Precision Manufacture Co., Ltd
No. 3001 Baohe Road, Baolong Industrial City, Longgang, Shenzhen, China .
Address: Tel: 0755-89888888
Tel: Fax:
Fax: EC Rep: Wellkang Ltd
16 Castle Street, Dover, CT16 1PW, England, UK

We, the manufacturer, declare under our sole responsibility that

the medical device(s)	Product Name	Single-use face mask
	Type/model, identification of product allowing traceability (Where applicable)	Flat/17.5cm×9.5cm
of class	according to annex IX of directive 93/42/EEC	Class I Medical Device (non-sterile, without measuring function)

is/are in conformity with the relevant provisions and requirements of directive 93/42/EEC, as amended by Directive 2007/47/EC.

Applied harmonised standards, national standards or other normative documents

EN 14683:2005
EN ISO 15223-1:2016
EN 1041:2008
EN ISO 14971:2012
EN ISO 10993-1:2009
EN ISO 10993-5:2009

Conformity assessment procedure


Module A (EC Declaration of Conformity (Annex VII) + Technical Files)

Notified Body
(name & number)
Certificate & number

NOT applicable

Signed on: 9 March 2020. Place: Shenzhen, Guangdong, China

Signature (on behalf of the manufacturer):  Liu Xiaoliang

Name of authorized signatory:  Liu Xiaoliang

Position held in the company: General Manager

Official Seal: